

Dr. Saracoglu graduated from Medicine Faculty of Istanbul University in 1988. He started his career in an intensive care unit at the International Hospital, İstanbul. He specialized in dermatovenerology at the Vakıf Gureba Hospital in İstanbul. His main interests include dermatologic surgery, dermatologic laser procedures, skin surface microscopy and laser gynecology procedures. Dr. Saracoglu is the founder of Estethica Medical Center in İstanbul.



#### Clinical Bulletin

J. LA&HA, Vol. 2014, No. 1; p. B10.

## Treatment of Post-Inflammatory Hyperpigmentation using Q-Switched Nd:YAG

Serafettin Saracoglu

### Parameters:

	QS- Nd:YAG (1064 nm)			
	Treatment 1	Treatment 2	Treatment 3	Treatment 4
Fluence:	2.5 J/cm <sup>2</sup>	2.9 J/cm <sup>2</sup>	3.1 J/cm <sup>2</sup>	3.3 J/cm <sup>2</sup>
Spot size:	6 mm	6 mm	6 mm	6 mm
Frequency:	8 Hz	8 Hz	8 Hz	8 Hz

### Treatment procedure:

A 32-year-old male patient with post-inflammatory hyperpigmentation characterized by irregular brown patches with poorly defined border, on the left cheek, came to our clinic. Post inflammatory hyperpigmentation had developed due to a trauma.

To treat the post inflammatory hyperpigmentation, a low-fluence Q-Switched Nd:YAG laser was used. A series of four treatments in 2-week intervals was needed to significantly reduce the hyperpigmented lesion, starting with a low fluence of 2.5 J/cm<sup>2</sup> in the first treatment and increasing it with the following treatments according to the desired endpoint effect. After the final Nd:YAG treatment, an additional pass using 532 nm Q-S KTP with 3 mm spot size, fluence 1.8 J/cm<sup>2</sup> and frequency of 2 Hz was performed to reduce residual epidermal spots. No anesthesia was used during the treatment.



Published by the Laser and Health Academy. All rights reserved. © 2014

Disclaimer: The intent of this Laser and Health Academy publication is to facilitate an exchange of information on the views, research results, and clinical experiences within the medical laser community. The contents of this publication are the sole responsibility of the authors and may not in any circumstances be regarded as official product information by the medical equipment manufacturers. When in doubt please check with the manufacturers whether a specific product or application has been approved or cleared to be marketed and sold in your country.



Before



After 4 treatments